

Dynamic nuclear polarisation in GaAs/AlGaAs dots observed with a low temperature confocal microscope - the attoCFM I.

We study nuclear spin effects in individual quantum dots pumped by circularly polarised excitation. Spin is transferred from optically generated electrons to the system of about 10.000 nuclear spins. The effect of nuclear spin polarisation on the electron spin can be described in terms of effective nuclear magnetic field (also referred to as Overhauser field). We observe optically generated Overhauser fields of the order of several Tesla localised in a single quantum dot of $4 \times 20 \times 20 \text{ nm}^3$ in size (see Figure 1).

Optically pumped nuclear polarisation has striking dynamic properties: for example in electron-charged InP/GaInP dots a high nuclear spin polarisation can be excited in a few ms, but survives after the light has been switched off for hundreds of seconds with some background polarisation observed after two hours!

Obviously, such measurement requires exceptional set-up stability, since for each point in the decay curve (see figure 3) a single dot is excited and probed by an individual pump-probe pulse pair. To remove uncertainty introduced by the noise in the CCD detectors, the whole decay curve has to be measured several times in the same conditions. The experiment is sometimes run for 25 hours, during which time the sample drift is considerably less than the 400 nm size of the aperture in the metallic shadow film through which the dot is optically excited.

Additional measurements of nuclear spin effects in charged and neutral quantum dots have been carried out with this setup (see [2]). These experiments only became possible after the attoCFM I microscope system had been introduced in our lab in November 2007, due to the following advantages:

- free space access conserving polarisation of light;
- exceptional mechanical stability both at zero and high magnetic fields;
- very small sample drift when scanning the field ($\ll 400 \text{ nm}$ up to 2.5 T).

The data has been generously provided by E. Chekhovich, M. Makhonin and A. Tartakovskii from the Department of Physics, Sheffield University, UK.

[1] D. Gammon, et al., Science 277, 85 (1997)

[2] E. Chekhovich, et al., submitted to PRL

RELATED PRODUCTS	
attoCFM I	highly flexible, low temperature confocal microscope
ANPxyz101/LT	high precision, piezo electric, inertial positioner
ANSxy100/LT	low temperature compatible xy-scanner
ASC500	SPM controller




Fig. 1: Nuclear spin polarisation in interface GaAs quantum dots. Very detailed power scans (1-2 hour scans) have been taken.




Fig. 2: Optically detected NMR (ODNMR) spectra. This allows insight into the local material composition on the nano-scale. This has been demonstrated only once before in 1997 by Gammon and co-workers [1]. The quality of data is considerably higher thanks to the stability of the attoCFM I.




Fig. 3: Nuclear polarisation decay measurements in GaAs dots. A very smooth decay curve is observed, with the nuclear spin decay most probably caused by the nuclear spin diffusion in the semiconductor matrix surrounding the dot.